

Datasheet

Dityrosine monoclonal antibody, clone 7D4 (FITC)

Catalog Number: MAB17421

Regulation Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody raised against synthetic Dityrosine (DT).

Clone Name: 7D4

Immunogen: Synthetic Dityrosine conjugated to Keyhole Limpet Hemocyanin (KLH).

Host: Mouse

Applications: ELISA, ICC, IF, WB
(See our web site product page for detailed applications information)

Protocols: See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Form: Liquid

Conjugation: FITC

Purification: Protein G purification

Isotype: IgG1

Recommend Usage: ELISA (1:1000)

Immunocytochemistry (1:50)

Immunofluorescence (1:50)

Western Blot (1:1000)

The optimal working dilution should be determined by the end user.

Storage Buffer: In PBS, pH 7.4 (50% glycerol, 0.09% sodium azide).

Storage Instruction: Store at -20°C.