

Datasheet

KDELRL1 monoclonal antibody, clone KR-10 (PerCP)

Entrez GeneID: 618184

Gene Symbol: KDELRL1

Catalog Number: MAB18056

Gene Alias: MGC139648

Regulation Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody raised against synthetic peptide of bovine KDELRL1.

Clone Name: KR-10

Immunogen: A synthetic peptide (conjugated with KLH) corresponding to amino acids 192-212 of bovine KDELRL1.

Host: Mouse

Reactivity: Bovine

Applications: ICC, IF, IHC-P, IP, WB
(See our web site product page for detailed applications information)

Protocols: See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Form: Liquid

Conjugation: PerCP

Purification: Protein G purification

Isotype: IgG1

Recommend Usage: Immunocytochemistry (1:1000)

Immunofluorescence (1:1000)

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections)

Immunoprecipitation

Western Blot (1:1000)

The optimal working dilution should be determined by the end user.

Storage Buffer: In PBS, pH 7.2 (50% glycerol, 0.09% sodium azide).

Storage Instruction: Store at -20°C.