

## Datasheet

### Acetylated Lysine monoclonal antibody, clone 7F8 (PerCP)

**Catalog Number:** MAB18080

**Regulation Status:** For research use only (RUO)

**Product Description:** Mouse monoclonal antibody raised against Acetylated Lysine.

**Clone Name:** 7F8

**Immunogen:** Acetylated KLH.

**Host:** Mouse

**Reactivity:** Mouse

**Applications:** ELISA, ICC, IF, IHC-P, IP, WB  
(See our web site product page for detailed applications information)

**Protocols:** See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

**Form:** Liquid

**Conjugation:** PerCP

**Purification:** Protein G purification

**Isotype:** IgG1

**Recommend Usage:** ELISA

Immunocytochemistry

Immunofluorescence

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) (1:100)

Immunoprecipitation

Western Blot (1:1000)

The optimal working dilution should be determined by the end user.

**Storage Buffer:** In PBS, pH 7.4 (50% glycerol, 0.09% sodium azide).

**Storage Instruction:** Store at -20°C.