

Datasheet

Nitrotyrosine monoclonal antibody, clone 39B6 (PerCP)

Storage Instruction: Store at -20°C.

Catalog Number: MAB18092

Regulation Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody raised against Nitrotyrosine.

Clone Name: 39B6

Immunogen: Hybridoma line 39B6.

Host: Mouse

Reactivity: Human, Mouse, Rat

Applications: AM, ELISA, Flow Cyt, ICC, IF, IHC-P, IP, WB

(See our web site product page for detailed applications information)

Protocols: See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Form: Liquid

Conjugation: PerCP

Purification: Protein G purification

Isotype: IgG2a

Recommend Usage: Antibody Microarray

ELISA

Flow Cytometry

Immunocytochemistry

Immunofluorescence

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) (1:100)

Immunoprecipitation

Western Blot (1:1400)

The optimal working dilution should be determined by the end user.

Storage Buffer: In PBS (50% glycerol, 0.09% sodium azide).