

## Datasheet

### Mmp9 monoclonal antibody, clone S51-82 (PerCP)

**Entrez GeneID:** 81687

**Gene Symbol:** Mmp9

**Catalog Number:** MAB18209

**Gene Alias:** -

**Regulation Status:** For research use only (RUO)

**Product Description:** Mouse monoclonal antibody raised against full length recombinant rat Mmp9.

**Clone Name:** S51-82

**Immunogen:** Recombinant protein corresponding to full length rat Mmp9.

**Host:** Mouse

**Reactivity:** Mouse, Rat

**Applications:** Flow Cyt, ICC, IF, IHC-P, IP, WB-Ti  
(See our web site product page for detailed applications information)

**Protocols:** See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

**Form:** Liquid

**Conjugation:** PerCP

**Purification:** Protein G purification

**Isotype:** IgG2a

**Recommend Usage:** Flow Cytometry  
Immunocytochemistry (1:100)  
Immunofluorescence (1:100)  
Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections)  
Immunoprecipitation  
Western Blot (1:1000)

The optimal working dilution should be determined by the end user.

**Storage Buffer:** In PBS, pH 7.4 (50% glycerol, 0.09% sodium azide).

**Storage Instruction:** Store at -20°C.