

Datasheet

Slc17a8 monoclonal antibody, clone S34-34 (PerCP)

Gene Symbol: Slc17a8

Gene Alias: Vglut3

Catalog Number: MAB18240

Regulation Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody raised against partial recombinant rat Slc17a8.

Clone Name: S34-34

Immunogen: Recombinant protein corresponding to amino acids 546-588 at C-terminus of rat Slc17a8.

Host: Mouse

Reactivity: Rat

Applications: ICC, IF, IHC-P, WB-Ti
(See our web site product page for detailed applications information)

Protocols: See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Form: Liquid

Conjugation: PerCP

Purification: Protein G Purification

Isotype: IgG1

Recommend Usage: Immunocytochemistry (1:100)
Immunofluorescence (1:100)
Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections)
Western Blot (1:1000)
The optimal working dilution should be determined by the end user.

Storage Buffer: In PBS, pH 7.4 (50% glycerol, 0.09% sodium azide).

Storage Instruction: Store at -20°C.

Entrez GeneID: 266767