

Datasheet

Gabbr1 monoclonal antibody, clone S93A-49 (PerCP)

Gene Symbol: Gabbr1

Gene Alias: -

Catalog Number: MAB18285

Regulation Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody raised against partial recombinant rat Gabbr1.

Clone Name: S93A-49

Immunogen: Recombinant protein corresponding to amino acids 873-977 at C-terminus of rat Gabbr1.

Host: Mouse

Reactivity: Rat

Applications: ICC, IF, WB-Ti

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Form: Liquid

Conjugation: PerCP

Purification: Protein G Purification

Isotype: IgG1

Recommend Usage: Immunocytochemistry

Immunofluorescence

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections)

Western Blot (1:1000)

The optimal working dilution should be determined by the end user.

Storage Buffer: In PBS, pH 7.4 (50% glycerol, 0.09% sodium azide).

Storage Instruction: Store at -20°C.

Entrez GeneID: 81657