

Datasheet

SLC1A1 monoclonal antibody, clone S180-41 (PerCP)

Gene Alias: EAAC1, EAAT3

Catalog Number: MAB18296

Regulation Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody raised against full length recombinant human SLC1A1.

Clone Name: S180-41

Immunogen: Recombinant protein corresponding to full length human SLC1A1.

Host: Mouse

Reactivity: Human, Rat

Applications: ICC, IF, IP, WB-Ti

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Form: Liquid

Conjugation: PerCP

Purification: Protein G Purification

Isotype: IgG1

Recommend Usage: Immunocytochemistry (1:100)

Immunofluorescence (1:100)

Immunoprecipitation

Western Blot (1:1000)

The optimal working dilution should be determined by the end user.

Storage Buffer: In PBS, pH 7.4 (50% glycerol, 0.09% sodium azide).

Storage Instruction: Store at -20°C.

Entrez GeneID: 6505

Gene Symbol: SLC1A1