

## Datasheet

### Grm1/Grm5 monoclonal antibody, clone S75-33 (PerCP)

**Gene Symbol:** Grm1

**Gene Alias:** Gprc1a

**Catalog Number:** MAB18308

**Regulation Status:** For research use only (RUO)

**Product Description:** Mouse monoclonal antibody raised against partial recombinant rat Grm1/Grm5.

**Clone Name:** S75-33

**Immunogen:** Recombinant protein corresponding to amino acids 824-1203 at C-terminus of rat Grm5

**Host:** Mouse

**Reactivity:** Rat

**Applications:** ICC, IF, IHC-P, IP, WB-Ti  
(See our web site product page for detailed applications information)

**Protocols:** See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

**Form:** Liquid

**Conjugation:** PerCP

**Purification:** Protein G Purification

**Isotype:** IgG1

**Recommend Usage:** Immunocytochemistry (1:100)  
Immunofluorescence (1:100)  
Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) (1:1000)  
Immunoprecipitation  
Western Blot (1:1000)

The optimal working dilution should be determined by the end user.

**Storage Buffer:** In PBS, pH 7.4 (50% glycerol, 0.09% sodium azide).

**Storage Instruction:** Store at -20°C.

**Entrez GeneID:** 24414|24418