

Datasheet

Abcc8 monoclonal antibody, clone S289-16 (PerCP)

Gene Symbol: Abcc8

Gene Alias: Sur, Sur1

Catalog Number: MAB18332

Regulation Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody raised against partial recombinant rat Abcc8.

Clone Name: S289-16

Immunogen: Recombinant protein corresponding to amino acids 1548-1582 at C-terminus of rat Abcc8.

Host: Mouse

Reactivity: Rat

Applications: ICC, IF, IHC-P, WB-Ti
(See our web site product page for detailed applications information)

Protocols: See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Form: Liquid

Conjugation: PerCP

Purification: Protein G Purification

Isotype: IgG1

Recommend Usage: Immunocytochemistry (1:100)
Immunofluorescence (1:100)
Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) (1:1000)
Western Blot (1:1000)

The optimal working dilution should be determined by the end user.

Storage Buffer: In PBS, pH 7.4 (50% glycerol, 0.09% sodium azide).

Storage Instruction: Store at -20°C.

Entrez GeneID: 25559