

## Datasheet

### Git1 monoclonal antibody, clone S39B-8 (PerCP)

Gene Alias: -

**Catalog Number:** MAB18356

**Regulation Status:** For research use only (RUO)

**Product Description:** Mouse monoclonal antibody raised against partial recombinant rat Git1.

**Clone Name:** S39B-8

**Immunogen:** Recombinant protein corresponding to amino acids 375-770 at C-terminus of rat Git1.

**Host:** Mouse

**Reactivity:** Rat

**Applications:** ICC, IF, IP, WB-Ti

(See our web site product page for detailed applications information)

**Protocols:** See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

**Form:** Liquid

**Conjugation:** PerCP

**Purification:** Protein G Purification

**Isotype:** IgG1

**Recommend Usage:** Immunocytochemistry

Immunofluorescence

Immunoprecipitation

Western Blot (1:1000)

The optimal working dilution should be determined by the end user.

**Storage Buffer:** In PBS, pH 7.4 (50% glycerol, 0.09% sodium azide).

**Storage Instruction:** Store at -20°C.

**Entrez GeneID:** 83709

**Gene Symbol:** Git1