

Datasheet

Nalcn monoclonal antibody, clone S187-7 (PerCP)

Gene Alias: Vgcnl1

Catalog Number: MAB18380

Regulation Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody raised against partial recombinant rat Nalcn.

Clone Name: S187-7

Immunogen: Recombinant protein corresponding to amino acids 1659-1738 at C-terminus of rat Nalcn.

Host: Mouse

Reactivity: Rat

Applications: ICC, IF, WB-Ti

(See our web site product page for detailed applications information)

Protocols: See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Form: Liquid

Conjugation: PerCP

Purification: Protein G Purification

Isotype: IgG1

Recommend Usage: Immunocytochemistry (1:100)

Immunofluorescence (1:100)

Western Blot (1:1000)

The optimal working dilution should be determined by the end user.

Storage Buffer: In PBS, pH 7.4 (50% glycerol, 0.09% sodium azide).

Storage Instruction: Store at -20°C.

Entrez GeneID: 266760

Gene Symbol: Nalcn