

Datasheet

BrdU monoclonal antibody, clone 131-14871 (AH4H7-1)

Catalog Number: MAB3900

Regulatory Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody raised against BrdU.

Clone Name: 131-14871 (AH4H7-1)

Immunogen: Bromodeoxyuridine (BrdU).

Host: Mouse

Applications: Flow Cyt, IHC
(See our web site product page for detailed applications information)

Protocols: See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Specificity: Bromodeoxyuridine (BrdU), a thymidine analogue, is incorporated into DNA during replication. Anti-BrdU reacts with cells that incorporated BrdU during proliferation activity in a variety of specimens such as tumors, blood, cell lines and tissues.

Form: Liquid

Isotype: IgG1, kappa

Recommend Usage: The optimal working dilution should be determined by the end user.

Storage Buffer: In 10 mM phosphate buffer, 150 mM NaCl, pH 7.4 (0.09% sodium azide)

Storage Instruction: Store at 4°C. For long term storage store at -20°C.
Aliquot to avoid repeated freezing and thawing.