

## Datasheet

### **c-Myc tag monoclonal antibody, clone 289-19510**

**Catalog Number:** MAB4164

**Regulatory Status:** For research use only (RUO)

**Product Description:** Mouse monoclonal antibody  
raised against synthetic peptides of c-Myc tag.

**Clone Name:** 289-19510

**Immunogen:** A synthetic peptide corresponding to  
c-Myc tag.

**Host:** Mouse

**Applications:** IHC-Fr, WB  
(See our web site product page for detailed applications  
information)

**Protocols:** See our web site at  
<http://www.abnova.com/support/protocols.asp> or product  
page for detailed protocols

**Specificity:** Myc tag (C-myc).

**Form:** Liquid

**Isotype:** IgG1, kappa

**Recommend Usage:** The optimal working dilution  
should be determined by the end user.

**Storage Buffer:** In 10 mM phosphate buffer, 150 mM  
NaCl, pH 7.4 (0.09% sodium azide)

**Storage Instruction:** Store at 4°C. For long term  
storage store at -80°C.  
Aliquot to avoid repeated freezing and thawing.