

Datasheet

CD4 monoclonal antibody, clone MEM-241 (PerCP)

Catalog Number: MAB4514

Regulatory Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody raised against partial recombinant CD4.

Clone Name: MEM-241

Immunogen: Recombinant Fc fusion protein corresponding to partial human CD4.

Host: Mouse

Theoretical MW (kDa): 55

Reactivity: Human

Applications: Flow Cyt, WB
(See our web site product page for detailed applications information)

Protocols: See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Specificity: This antibody recognizes CD4 antigen, a 55 KDa transmembrane glycoprotein expressed on a subset of T lymphocytes (helper T-cells) and also on monocytes, tissue macrophages and granulocytes.

Form: Liquid

Conjugation: PerCP

Isotype: IgG1

Recommend Usage: Flow Cytometry (10 ug/mL)
The optimal working dilution should be determined by the end user.

Storage Buffer: In PBS (0.2% BSA, 0.09% sodium azide)

Storage Instruction: Store in the dark at 4°C. Do not freeze.
Avoid prolonged exposure to light.

Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 920

Gene Symbol: CD4

Gene Alias: CD4mut

References:

1. CD molecules 2006--human cell differentiation molecules. Zola H, Swart B, Banham A, Barry S, Beare A, Bensussan A, Boumsell L, D Buckley C, Buhring HJ, Clark G, Engel P, Fox D, Jin BQ, Macardle PJ, Malavasi F, Mason D, Stockinger H, Yang X. J Immunol Methods. 2007 Jan 30;319(1-2):1-5. Epub 2006 Dec 4.
2. LIME: a new membrane Raft-associated adaptor protein involved in CD4 and CD8 coreceptor signaling. Brdickova N, Brdicka T, Angelisova P, Horvath O, Spicka J, Hilgert I, Paces J, Simeoni L, Kliche S, Merten C, Schraven B, Horejsi V. J Exp Med. 2003 Nov 17;198(10):1453-62. Epub 2003 Nov 10.
3. Cell surface receptors, virus entry and tropism of primate lentiviruses. Clapham PR, McKnight A. J Gen Virol. 2002 Aug;83(Pt 8):1809-29.