

Datasheet

Ptprc monoclonal antibody, clone ANK74

Catalog Number: MAB5278

Regulatory Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody raised against native Ptprc.

Clone Name: ANK74

Immunogen: Native purified Ptprc from IL-2-activated cultured Wag rats NK cells.

Host: Mouse

Reactivity: Rat

Applications: Flow Cyt, IHC-Fr
(See our web site product page for detailed applications information)

Protocols: See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Specificity: This is a rat strain-independent antibody to CD45. This monoclonal antibody recognizes a common epitope of rat CD45, present on all rat leucocytes.

Form: Liquid

Isotype: IgG1

Recommend Usage: The optimal working dilution should be determined by the end user.

Storage Buffer: In culture supernatant (0.05% sodium azide)

Storage Instruction: Store at 4°C. For long term storage store at -20°C.
Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 24699

Gene Symbol: Ptprc

Gene Alias: CD45, L-CA, Lca, RT7, T200