

Datasheet

Influenza A virus nucleoprotein monoclonal antibody, clone 431 (FITC)

Catalog Number: MAB5468

Regulatory Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody
raised against influenza A virus nucleoprotein

Clone Name: 431

Immunogen: Native purified Influenza A virus virions.

Host: Mouse

Reactivity: Viruses

Applications: IF, IHC

(See our web site product page for detailed applications
information)

Protocols: See our web site at
<http://www.abnova.com/support/protocols.asp> or product
page for detailed protocols

Specificity: Influenza A-specific for nucleoprotein.
Reactive with all A strains.

Form: Liquid

Conjugation: FITC

Isotype: IgG2a

Recommend Usage: Immunofluorescence (1:5-1:20)
Immunohistochemistry (1:5-1:20)
The optimal working dilution should be determined by
the end user.

Storage Buffer: In 10 mM PBS, pH 7.2 (10 mg/mL BSA,
0.09% sodium azide)

Storage Instruction: Store at -20°C.
Aliquot to avoid repeated freezing and thawing.