

Datasheet

Human immunodeficiency virus type 1 p24 monoclonal antibody, clone 491 (FITC)

Catalog Number: MAB5472

Regulatory Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody raised against native human immunodeficiency virus type 1 p24.

Clone Name: 491

Immunogen: Native purified human immunodeficiency virus type 1 (HIV-1) p24.

Host: Mouse

Reactivity: Viruses

Applications: EIA, ELISA, IF, IHC, WB
(See our web site product page for detailed applications information)

Protocols: See our web site at
<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Specificity: Specific for HIV1 p24. Does not cross-react with other HIV polypeptides.

Form: Liquid

Conjugation: FITC

Isotype: IgG1

Recommend Usage: The optimal working dilution should be determined by the end user.

Storage Buffer: In 10 mM PBS, pH 7.2 (10 mg/mL BSA, 0.09% sodium azide)

Storage Instruction: Store at -20°C.
Aliquot to avoid repeated freezing and thawing.