

## Datasheet

### Respiratory syncytial virus Nucleoprotein monoclonal antibody, clone 671 (FITC)

**Catalog Number:** MAB5510

**Regulatory Status:** For research use only (RUO)

**Product Description:** Mouse monoclonal antibody raised against Respiratory syncytial virus Nucleoprotein.

**Clone Name:** 671

**Immunogen:** Native from respiratory syncytial virus viral lysate (Long strain).

**Host:** Mouse

**Reactivity:** Viruses

**Applications:** ELISA, IF  
(See our web site product page for detailed applications information)

**Protocols:** See our web site at  
<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

**Specificity:** Specific for the nucleoprotein of RSV.

**Form:** Liquid

**Conjugation:** FITC

**Isotype:** IgG2a

**Recommend Usage:** The optimal working dilution should be determined by the end user.

**Storage Buffer:** In 10 mM PBS, pH 7.2 (10 mg/mL BSA, 0.09% sodium azide)

**Storage Instruction:** Store at -20°C.  
Aliquot to avoid repeated freezing and thawing.