

Datasheet

MYC monoclonal antibody, clone Myc.A7

Catalog Number: MAB5543

Regulatory Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody raised against synthetic peptide of MYC.

Clone Name: Myc.A7

Immunogen: A synthetic peptide (conjugated with KLH) corresponding to amino acids 410-419 of human MYC.

Sequence: EQKLISEEDL

Host: Mouse

Reactivity: Human

Applications: IHC, IP, WB

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Specificity: This antibody recognizes overexpressed proteins containing MYC epitope tag fused to either amino- or carboxy-termini of targeted proteins in transfected mammalian cells.

Form: Liquid

Isotype: IgG1

Recommend Usage: Western Blot (1:1000)

immunoprecipitation (1:200)

Immunohistochemistry (1:200)

The optimal working dilution should be determined by the end user.

Storage Buffer: In PBS, pH 7.4 (0.05% sodium azide)

Storage Instruction: Store at 4°C. For long term storage store at -20°C.

Aliquot to avoid repeated freezing and thawing.