

Datasheet

Cd19 monoclonal antibody, clone 6D5 (PE/Cy5.5)

Catalog Number: MAB5693

Regulatory Status: For research use only (RUO)

Product Description: Rat monoclonal antibody raised against Cd19.

Clone Name: 6D5

Immunogen: Mouse CD19-expressing K562 human erythroleukemia cells.

Host: Rat

Reactivity: Mouse

Applications: Flow Cyt, IP
(See our web site product page for detailed applications information)

Protocols: See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Specificity: CD19 B-cell differentiation antigen, Mr 95 KDa.

Form: Liquid

Conjugation: PE/Cy5.5

Isotype: IgG2a, kappa

Recommend Usage: Flow Cytometry (0.2 ug/10⁶ cells)
The optimal working dilution should be determined by the end user.

Storage Buffer: In PBS (0.09% sodium azide)

Storage Instruction: Store in the dark at 4°C. Do not freeze.

Avoid prolonged exposure to light.

Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 12478

Gene Symbol: Cd19

Gene Alias: AW495831

References:

1. The signaling activity of murine CD19 is regulated during cell development. Krop I, Shaffer AL, Fearon DT, Schlissel MS. J Immunol. 1996 Jul 1;157(1):48-56.
2. Self-renewal of B-1 lymphocytes is dependent on CD19. Krop I, de Fougerolles AR, Hardy RR, Allison M, Schlissel MS, Fearon DT. Eur J Immunol. 1996 Jan;26(1):238-42.
3. The CD19-CR2-TAPA-1 complex, CD45 and signaling by the antigen receptor of B lymphocytes. Fearon DT. Curr Opin Immunol. 1993 Jun;5(3):341-8.