

Datasheet

I-Ek monoclonal antibody, clone 17-3-3S (FITC)

DH. J Immunol. 1980 Feb;124(2):533-40.

Catalog Number: MAB5804**Regulatory Status:** For research use only (RUO)**Product Description:** Mouse monoclonal antibody raised against I-Ek.**Clone Name:** 17-3-3S**Immunogen:** Native purified I-Ek from C3H mouse skin graft and splenocytes.**Host:** Mouse**Reactivity:** Mouse**Applications:** Flow Cyt
(See our web site product page for detailed applications information)**Protocols:** See our web site at
<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols**Specificity:** An epitope on mouse I-Ek MHC class II alloantigen.**Form:** Liquid**Conjugation:** FITC**Isotype:** IgG2a, kappa**Recommend Usage:** Flow Cytometry (1 ug/10⁶ cells)
The optimal working dilution should be determined by the end user.**Storage Buffer:** In PBS (0.09% sodium azide)**Storage Instruction:** Store in the dark at 4°C. Do not freeze.

Avoid prolonged exposure to light.

Aliquot to avoid repeated freezing and thawing.

References:

1. Hybridoma cell lines secreting monoclonal antibodies to mouse H-2 and Ia antigens. Ozato K, Mayer N, Sachs