

Datasheet

H-2Kk monoclonal antibody, clone 16-3-22S (PE)

Catalog Number: MAB5831

Regulatory Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody raised against native H-2Kk.

Clone Name: 16-3-22S

Immunogen: Native purified from C3H.SW mouse splenocytes.

Host: Mouse

Reactivity: Mouse

Applications: Flow Cyt, IHC-Fr
(See our web site product page for detailed applications information)

Protocols: See our web site at
<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Specificity: An epitope on mouse H-2Kk MHC class I alloantigen.

Form: Liquid

Conjugation: PE

Isotype: IgG2a, kappa

Recommend Usage: Flow Cytometry (0.2 ug/10⁶ cells)
The optimal working dilution should be determined by the end user.

Storage Buffer: In PBS (0.09% sodium azide)

Storage Instruction: Store in the dark at 4°C. Do not freeze.

Avoid prolonged exposure to light.

Aliquot to avoid repeated freezing and thawing.