

Datasheet

H-2Kk monoclonal antibody, clone 36-7-5 (FITC)

Catalog Number: MAB5832

Regulatory Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody raised against native H-2Kk.

Clone Name: 36-7-5

Immunogen: Native purified from A.AL mouse splenocytes.

Host: Mouse

Reactivity: Mouse

Applications: Flow Cyt, IHC-Fr
(See our web site product page for detailed applications information)

Protocols: See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Specificity: An epitope on mouse H-2Kk MHC class I alloantigen.

Form: Liquid

Conjugation: FITC

Isotype: IgG2a, kappa

Recommend Usage: Flow Cytometry (1 ug/10⁶ cells)
The optimal working dilution should be determined by the end user.

Storage Buffer: In PBS (0.09% sodium azide)

Storage Instruction: Store in the dark at 4°C. Do not freeze.

Avoid prolonged exposure to light.

Aliquot to avoid repeated freezing and thawing.