

Datasheet

CD8A/CD8B monoclonal antibody, clone RFT-8 (PE/TXRD)

Catalog Number: MAB6024

Regulatory Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody raised against native CD8A/CD8B.

Clone Name: RFT-8

Immunogen: Native purified from human thymocytes and E-rosetted lymphocytes.

Host: Mouse

Reactivity: Human

Applications: Flow Cyt, IHC-Fr
(See our web site product page for detailed applications information)

Protocols: See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Specificity: Specificity human CD8

Form: Liquid

Conjugation: PE/TXRD

Isotype: IgG1

Recommend Usage: Flow Cytometry (10 ul/10⁶ cells)
The optimal working dilution should be determined by the end user.

Storage Buffer: In PBS (0.09% sodium azide)

Storage Instruction: Store in the dark at 4°C. Do not freeze.

Avoid prolonged exposure to light.

Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 925|926

Gene Symbol: CD8A

Gene Alias: CD8, Leu2, MAL, p32

References:

1. The CD8 coreceptor revisited: one chain good, two chains better. Zamoyska R. Immunity. 1994 Jul;1(4):243-6.