

## Datasheet

### GSK3A monoclonal antibody, clone SBGSK3a4 (SEPH)

**Catalog Number:** MAB6088

**Regulatory Status:** For research use only (RUO)

**Product Description:** Mouse monoclonal antibody raised against recombinant GSK3A.

**Clone Name:** SBGSK3a4

**Immunogen:** Recombinant protein corresponding to GSK3A.

**Host:** Mouse

**Reactivity:** Human, Mouse, Rat

**Applications:** IP, WB

(See our web site product page for detailed applications information)

**Protocols:** See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

**Specificity:** Specificity human, mouse and rat GSK-3alpha

**Form:** Liquid

**Conjugation:** PE/TXRD

**Isotype:** IgG2a

**Recommend Usage:** Immunoprecipitation (40 ul)

The optimal working dilution should be determined by the end user.

**Storage Buffer:** In 100 mM BBS, pH 8.2

**Storage Instruction:** Store in the dark at 4°C. Do not freeze.

Avoid prolonged exposure to light.

Aliquot to avoid repeated freezing and thawing.