

## Datasheet

### COL4 monoclonal antibody, clone 2F11 (FITC)

**Catalog Number:** MAB6091

**Regulatory Status:** For research use only (RUO)

**Product Description:** Mouse monoclonal antibody raised against native COL4.

**Clone Name:** 2F11

**Immunogen:** Native human COL4.

**Host:** Mouse

**Reactivity:** Human

**Applications:** ELISA, IHC-Fr  
(See our web site product page for detailed applications information)

**Protocols:** See our web site at  
<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

**Specificity:** Specificity Reacts with human Type IV collagen. Cross reactivity with other species has not been established.

**Form:** Liquid

**Conjugation:** FITC

**Isotype:** IgG1

**Recommend Usage:** Immunohistochemistry (2 ug/mL)  
The optimal working dilution should be determined by the end user.

**Storage Buffer:** In PBS (0.09% sodium azide)

**Storage Instruction:** Store in the dark at 4°C. Do not freeze.

Avoid prolonged exposure to light.

Aliquot to avoid repeated freezing and thawing.