

Datasheet

COL4 monoclonal antibody, clone 2F11 (PE)

Catalog Number: MAB6093

Regulatory Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody raised against native COL4.

Clone Name: 2F11

Immunogen: Native human COL4.

Host: Mouse

Reactivity: Human

Applications: ELISA, IHC-Fr
(See our web site product page for detailed applications information)

Protocols: See our web site at
<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Specificity: Specificity Reacts with human Type IV collagen. Cross reactivity with other species has not been established.

Form: Liquid

Conjugation: PE

Isotype: IgG1

Recommend Usage: The optimal working dilution should be determined by the end user.

Storage Buffer: In PBS (0.09% sodium azide)

Storage Instruction: Store at 4°C. Do not freeze.