

Datasheet

Klf4 monoclonal antibody

Catalog Number: MAB8725

Regulatory Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody raised against full length recombinant Klf4.

Immunogen: Recombinant protein corresponding to full length mouse Klf4.

Host: Mouse

Reactivity: Human, Mouse

Applications: ELISA, WB-Ce
(See our web site product page for detailed applications information)

Protocols: See our web site at
<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Specificity: This antibody is specific to Klf4.

Form: Liquid

Purification: Protein G purification

Concentration: 1 mg/mL

Recommend Usage: Western Blot (1:100-1:2000)
ELISA (1:5000-1:20000)
The optimal working dilution should be determined by the end user.

Storage Buffer: In PBS, pH 7.4 (0.02% sodium azide, 1% BSA, 50% glycerol)

Storage Instruction: Store at 4°C for three months. For long term storage store at -20°C.
Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 16600

Gene Symbol: Klf4

Gene Alias: EZF, Gklf, Zie