

Datasheet

Sox2 monoclonal antibody

Catalog Number: MAB8752

Regulatory Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody raised against full length recombinant Sox2.

Immunogen: Recombinant protein corresponding to full length mouse Sox2.

Host: Mouse

Reactivity: Mouse

Applications: ELISA, WB-Re
(See our web site product page for detailed applications information)

Protocols: See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Specificity: This antibody is specific to Sox2.

Form: Liquid

Purification: Protein G purification

Concentration: 1 mg/mL

Recommend Usage: Western Blot (1:100-1:2000)
ELISA (1:5000-1:20000)
The optimal working dilution should be determined by the end user.

Storage Buffer: In PBS, pH 7.4 (0.02% sodium azide, 1% BSA, 50% glycerol)

Storage Instruction: Store at 4°C for three months. For long term storage store at -20°C.
Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 20674

Gene Symbol: Sox2

Gene Alias: Sox-2, lcc, ysb