

Datasheet

POU5F1 monoclonal antibody, clone 97q8

Gene Symbol: POU5F1

Gene Alias: MGC22487, OCT3, OCT4, OTF3, OTF4

Catalog Number: MAB8811

Regulatory Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody raised against partial recombinant POU5F1.

Clone Name: 97q8

Immunogen: Recombinant protein corresponding to amino acids 102-265 of human POU5F1.

Host: Mouse

Reactivity: Human

Applications: ELISA, WB-Re

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Specificity: This antibody recognizes recombinant protein POU5F1.

Form: Liquid

Purification: Affinity purification

Concentration: 1 ug/uL

Isotype: IgG1

Recommend Usage: Western Blot (1:1000)

The optimal working dilution should be determined by the end user.

Storage Buffer: In citrate-Tris-HCl, pH 7.0 (0.02% Proclin 300)

Storage Instruction: Store at -20°C.

Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 5460