

## Datasheet

### **SARS coronavirus N Protein monoclonal antibody, clone 38j4**

**Catalog Number:** MAB8899

**Regulatory Status:** For research use only (RUO)

**Product Description:** Mouse monoclonal antibody raised against full length recombinant *SARS coronavirus* N Protein.

**Clone Name:** 38j4

**Amount:** 100 uL

**Immunogen:** Recombinant protein corresponding to full length *SARS coronavirus* N Protein.

**Host:** Mouse

**Reactivity:** Viruses

**Applications:** ELISA, WB-Re  
(See our web site product page for detailed applications information)

**Protocols:** See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

**Specificity:** This antibody recognizes recombinant protein SARS N Protein.

**Form:** Liquid

**Purification:** Affinity purification

**Concentration:** 1 ug/uL

**Isotype:** IgG1

**Recommend Usage:** Western Blot (1:1000)  
The optimal working dilution should be determined by the end user.

**Storage Buffer:** In citrate-Tris-HCl, pH7.0 (0.02% Proclin 300)

**Storage Instruction:** Store at -20°C.  
Aliquot to avoid repeated freezing and thawing.