

Datasheet

Itgb2 monoclonal antibody, clone M18/2 (FITC)

Catalog Number: MAB9797

Regulatory Status: For research use only (RUO)

Product Description: Rat monoclonal antibody raised against native Itgb2.

Clone Name: M18/2

Immunogen: Native purified Itgb2 from mouse cytotoxic T cell glycoproteins.

Host: Rat

Theoretical MW (kDa): 95

Reactivity: Mouse

Applications: Flow Cyt
(See our web site product page for detailed applications information)

Protocols: See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Specificity: The rat monoclonal antibody M18/2 recognizes Itgb2 (Cd18), a 95 kDa type I transmembrane protein expressed on all leukocytes.

Form: Liquid

Conjugation: FITC

Concentration: 0.5 mg/mL

Isotype: IgG2a

Recommend Usage: Flow Cytometry (3 ug/mL)
The optimal working dilution should be determined by the end user.

Storage Buffer: In PBS, pH 7.4 (0.09% sodium azide)

Storage Instruction: Store in the dark at 4°C. Avoid prolonged exposure to light. Do not freeze.

Entrez GeneID: 16414

Gene Symbol: Itgb2

Gene Alias: 2E6, AI528527, Cd18, LAD, LCAMB, Lfa1, MF17