

## Datasheet

### IgG (Fc) monoclonal antibody, clone EM-07 (FITC)

**Catalog Number:** MAB9827

**Regulatory Status:** For research use only (RUO)

**Product Description:** Mouse monoclonal antibody raised against IgG (Fc).

**Clone Name:** EM-07

**Immunogen:** Fusion protein of human IgG Fc fragment.

**Host:** Mouse

**Reactivity:** Human

**Applications:** IF

(See our web site product page for detailed applications information)

**Protocols:** See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

**Form:** Liquid

**Conjugation:** FITC

**Concentration:** 0.1 mg/mL

**Isotype:** IgG1

**Recommend Usage:** Immunofluorescence (2-5 ug/mL)

The optimal working dilution should be determined by the end user.

**Storage Buffer:** In PBS, pH 7.4 (0.02% BSA, 0.09% sodium azide)

**Storage Instruction:** Store in the dark at 4°C. Avoid prolonged exposure to light. Do not freeze.