

## Datasheet

### Itgam monoclonal antibody, clone M1/70 (FITC)

**Catalog Number:** MAB9861

**Regulatory Status:** For research use only (RUO)

**Product Description:** Rat monoclonal antibody raised against Itgam.

**Clone Name:** M1/70

**Immunogen:** Native Itgam from mouse splenic T cells and splenocytes.

**Host:** Rat

**Reactivity:** Mouse

**Applications:** Flow Cyt  
(See our web site product page for detailed applications information)

**Protocols:** See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

**Form:** Liquid

**Conjugation:** FITC

**Isotype:** IgG2b, kappa

**Recommend Usage:** Flow Cytometry (? 1 ug/10<sup>6</sup> cells)  
The optimal working dilution should be determined by the end user.

**Storage Buffer:** In PBS (0.09% sodium azide)

**Storage Instruction:** Store at 4°C. Do not freeze.

**Entrez GeneID:** 16409

**Gene Symbol:** Itgam

**Gene Alias:** CD11b/CD18, CR3, CR3A, Cd11b, F730045J24Rik, Ly-40, MAC1, Mac-1, Mac-1a