

Datasheet

Itga5 monoclonal antibody, clone 5H10-27 (FITC)

Catalog Number: MAB9872

Regulatory Status: For research use only (RUO)

Product Description: Rat monoclonal antibody raised against Itga5.

Clone Name: 5H10-27 (MFR5)

Immunogen: Native Itga5 from mouse mast cell line MC/9.

Host: Rat

Reactivity: Mouse

Applications: Flow Cyt
(See our web site product page for detailed applications information)

Protocols: See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Form: Liquid

Conjugation: FITC

Isotype: IgG2a, kappa

Recommend Usage: Flow Cytometry (? 1 ug/10⁶ cells)
The optimal working dilution should be determined by the end user.

Storage Buffer: In PBS (0.09% sodium azide)

Storage Instruction: Store at 4°C. Do not freeze.

Entrez GeneID: 16402

Gene Symbol: Itga5

Gene Alias: Cd49e, Fnra