

Datasheet

NANOG polyclonal antibody

Catalog Number: PAB0375

Regulatory Status: For research use only (RUO)

Product Description: Rabbit polyclonal antibody raised against synthetic peptide of NANOG.

Immunogen: A synthetic peptide corresponding to amino acids 1-50 of human NANOG.

Host: Rabbit

Reactivity: Chimpanzee, Goat, Human, Monkey, Mouse

Applications: WB-Ce

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Form: Liquid

Recommend Usage: The optimal working dilution should be determined by the end user.

Storage Buffer: In PBS (0.05% BSA, 0.05% sodium azide)

Storage Instruction: Store at 4°C. For long term storage store at -20°C.
Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 79923

Gene Symbol: NANOG

Gene Alias: -