

Datasheet

Goat Anti-Chicken IgG (H&L) secondary antibody (Biotin)

1. Protein biotinylation. Bayer EA, Wilchek M. Methods Enzymol. 1990;184:138-60.

Catalog Number: PAB10502

Regulatory Status: For research use only (RUO)

Product Description: Goat anti-chicken IgG recognizes chicken IgG whole molecule. This secondary antibody was purified using antigen affinity chromatography. The antibody is conjugated with Biotin.

Immunogen: Chicken IgG whole molecule

Host: Goat

Reactivity: Chicken

Applications: ELISA, IF, IHC, WB
(See our web site product page for detailed applications information)

Protocols: See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Form: Lyophilized

Conjugation: Biotin

Purification: This product was prepared from monospecific antiserum by immunoaffinity chromatography, followed by solid phase adsorption(s) to remove any unwanted reactivities.

Conjugation Note: 10-20 BAC molecules per goat IgG molecule

Recommend Usage: The optimal working dilution should be determined by the end user.

Storage Buffer: Lyophilized from 20 mM KH_2PO_4 , 150 mM NaCl, pH 7.2 (0.01% sodium azide)

Storage Instruction: Store at 4°C on dry atmosphere. After reconstitution with deionized water, store at -20°C or lower.
Aliquot to avoid repeated freezing and thawing.

References: