

## Datasheet

### Rabbit Anti-Chicken IgG F(ab')<sub>2</sub> secondary antibody (Peroxidase)

Immunol Methods. 1981;47(2):129-44.

**Catalog Number:** PAB10512

**Regulatory Status:** For research use only (RUO)

**Product Description:** Rabbit anti-chicken IgG recognizes chicken IgG F(ab')<sub>2</sub> fragment. This secondary antibody was purified using antigen affinity chromatography. The antibody is conjugated with Peroxidase.

**Immunogen:** Chicken IgG F(ab')<sub>2</sub> fragment

**Host:** Rabbit

**Reactivity:** Chicken

**Applications:** Dot, EIA, ELISA, IHC, WB  
(See our web site product page for detailed applications information)

**Protocols:** See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

**Form:** Lyophilized

**Conjugation:** HRP

**Purification:** This product was prepared from monospecific antiserum by immunoaffinity chromatography, followed by solid phase adsorption(s) to remove any unwanted reactivities.

**Recommend Usage:** The optimal working dilution should be determined by the end user.

**Storage Buffer:** Lyophilized from 20 mM KH<sub>2</sub>PO<sub>4</sub>, 150 mM NaCl, pH 7.2 (0.01% gentamicin sulfate)

**Storage Instruction:** Store at 4°C on dry atmosphere. After reconstitution with 1 mL of deionized water, store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

**References:**

1. Immunohistochemistry with enzyme labeled antibodies: a brief review. Farr AG, Nakane PK. J