

Datasheet

Rabbit Anti-Donkey IgG (H&L) secondary antibody (Texas Red)

microfluorometric and fluorescence microscopic studies.
Titus JA, Haugland R, Sharrow SO, Segal DM. J
Immunol Methods. 1982;50(2):193-204.

Catalog Number: PAB10561

Regulatory Status: For research use only (RUO)

Product Description: Rabbit anti-donkey IgG recognizes donkey IgG whole molecule. This secondary antibody was purified using antigen affinity chromatography. The antibody is conjugated with Texas Red.

Immunogen: Donkey IgG whole molecule

Host: Rabbit

Reactivity: Donkey

Applications: Flow Cyt, IF
(See our web site product page for detailed applications information)

Protocols: See our web site at
<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Form: Lyophilized

Conjugation: Texas Red™

Purification: This product was prepared from monospecific antiserum by immunoaffinity chromatography, followed by solid phase adsorption(s) to remove any unwanted reactivities.

Recommend Usage: The optimal working dilution should be determined by the end user.

Storage Buffer: Lyophilized from 20 mM KH₂PO₄, 150 mM NaCl, pH 7.2 (0.01% sodium azide)

Storage Instruction: Store at 4 °C on dry atmosphere. After reconstitution with 1 mL of deionized water, store at -20 °C or lower. Aliquot to avoid repeated freezing and thawing.

References:

1. Texas Red, a hydrophilic, red-emitting fluorophore for use with fluorescein in dual parameter flow