

Datasheet

Rabbit Anti-Goat IgG (H&L) secondary antibody (Peroxidase)

Immunol Methods. 1981;47(2):129-44.

Catalog Number: PAB10576**Regulatory Status:** For research use only (RUO)**Product Description:** Rabbit anti-goat IgG recognizes goat IgG whole molecule. This secondary antibody was purified using antigen affinity chromatography. The antibody is conjugated with Peroxidase.**Immunogen:** Goat IgG whole molecule**Host:** Rabbit**Reactivity:** Goat**Applications:** Dot, EIA, ELISA, IHC, WB
(See our web site product page for detailed applications information)**Protocols:** See our web site at
<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols**Form:** Lyophilized**Conjugation:** HRP**Purification:** This product was prepared from monospecific antiserum by immunoaffinity chromatography, followed by solid phase adsorption(s) to remove any unwanted reactivities.**Recommend Usage:** The optimal working dilution should be determined by the end user.**Storage Buffer:** Lyophilized from 20 mM KH₂PO₄, 150 mM NaCl, pH 7.2 (0.01% gentamicin) Lyophilized from sulfate)**Storage Instruction:** Store at 4°C on dry atmosphere.
After reconstitution with 1 mL of deionized water, store at -20°C or lower.
Aliquot to avoid repeated freezing and thawing.**References:**

1. Immunohistochemistry with enzyme labeled antibodies: a brief review. Farr AG, Nakane PK. J