

Datasheet

Rabbit Anti-Guinea pig IgG F(ab')₂ secondary antibody (Fluorescein)

Catalog Number: PAB10632

Regulatory Status: For research use only (RUO)

Product Description: Rabbit anti-guinea pig IgG recognizes Guinea pig IgG F(ab')₂ fragment. This secondary antibody was purified using antigen affinity chromatography. The antibody is conjugated with Fluorescein.

Immunogen: Guinea pig IgG F(ab')₂ fragment

Host: Rabbit

Reactivity: Guinea pig

Applications: Flow Cyt, IF
(See our web site product page for detailed applications information)

Protocols: See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Form: Lyophilized

Conjugation: FITC

Purification: This product was prepared from monospecific antiserum by immunoaffinity chromatography, followed by solid phase adsorption(s) to remove any unwanted reactivities.

Conjugation Note: 2.2 moles DTAF per mole of rabbit IgG

Recommend Usage: The optimal working dilution should be determined by the end user.

Storage Buffer: Lyophilized from 10 mM KH₂PO₄, 150 mM NaCl, pH 7.2 (0.01% thimerosal)

Storage Instruction: Store at 4°C on dry atmosphere. After reconstitution with 1 mL of deionized water, store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

References:

1. Protein biotinylation. Bayer EA, Wilchek M. Methods Enzymol. 1990;184:138-60.