

Datasheet

Rabbit Anti-Pig IgG (H&L) secondary antibody (Fluorescein)

Catalog Number: PAB10793

Regulatory Status: For research use only (RUO)

Product Description: Rabbit anti-pig IgG recognizes pig IgG whole molecule. This secondary antibody was purified using antigen affinity chromatography. The antibody is conjugated with Fluorescein.

Immunogen: Pig IgG whole molecule

Host: Rabbit

Reactivity: Pig

Applications: Flow Cyt, IF
(See our web site product page for detailed applications information)

Protocols: See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Form: Lyophilized

Conjugation: FITC

Purification: This product was prepared from monospecific antiserum by immunoaffinity chromatography, followed by solid phase adsorption(s) to remove any unwanted reactivities.

Conjugation Note: 2.4 moles DTAF per mole of rabbit IgG

Recommend Usage: The optimal working dilution should be determined by the end user.

Storage Buffer: Lyophilized from 5 mM sodium phosphate, 70 mM NaCl, 50% (v/v) glycerol, pH 7.4

Storage Instruction: Store at 4°C on dry atmosphere. After reconstitution with 1 mL of deionized water, store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

References:

1. Immunofluorescence using dichlorotriazinylaminofluorescein (DTAF). I. Preparation and fractionation of labelled IgG. Blakeslee D, Baines MG. J Immunol Methods. 1976;13(3-4):305-20.