

Datasheet

Goat Anti-Rabbit IgG F(ab')₂ secondary antibody (Fluorescein)

Catalog Number: PAB10827

Regulatory Status: For research use only (RUO)

Product Description: Goat anti-rabbit IgG recognizes rabbit IgG F(ab')₂ fragment. This secondary antibody was purified using antigen affinity chromatography. The antibody is conjugated with Fluorescein.

Immunogen: Rabbit IgG F(ab')₂ fragment

Host: Goat

Reactivity: Rabbit

Applications: Flow Cyt, IF
(See our web site product page for detailed applications information)

Protocols: See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Form: Lyophilized

Conjugation: FITC

Purification: This product was prepared from monospecific antiserum by immunoaffinity chromatography, followed by solid phase adsorption(s) to remove any unwanted reactivities.

Isotype: IgG

Conjugation Note: 5.3 moles Fluorescein (FITC) per mole of IgG

Recommend Usage: The optimal working dilution should be determined by the end user.

Storage Buffer: Lyophilized from 0.01 M Sodium Phosphate, 0.15 M Sodium Chloride, pH 7.2 (0.01% (w/v) thimerosal, 10 mg/mL BSA (Immunoglobulin and Protease free))

Storage Instruction: Store at 4°C on dry atmosphere. After reconstitution with deionized water, store at -20°C

or lower.

Aliquot to avoid repeated freezing and thawing.

References:

1. Immunofluorescence using dichlorotriazinylaminofluorescein (DTAF). I. Preparation and fractionation of labelled IgG. Blakeslee D, Baines MG. J Immunol Methods. 1976;13(3-4):305-20.