

Datasheet

Chicken Anti-Rat IgG (H&L) secondary antibody (Biotin)

Catalog Number: PAB10848

Regulatory Status: For research use only (RUO)

Product Description: Chicken anti-rat IgG recognizes rat IgG whole molecule. This secondary antibody was purified using antigen affinity chromatography. The antibody is conjugated with Biotin.

Immunogen: Rat IgG whole molecule

Host: Chicken

Reactivity: Rat

Applications: ELISA, IF, IHC, WB
(See our web site product page for detailed applications information)

Protocols: See our web site at
<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Form: Lyophilized

Conjugation: Biotin

Purification: This product was prepared from monospecific antiserum by immunoaffinity chromatography, followed by solid phase adsorption(s) to remove any unwanted reactivities.

Isotype: IgG

Conjugation Note:

Recommend Usage: The optimal working dilution should be determined by the end user.

Storage Buffer: Lyophilized from 20 mM potassium phosphate, 150 mM NaCl, pH 7.2 (50% glycerol, 0.01% sodium azide, 10 mg/mL BSA, Immunoglobulin and Protease free).

Storage Instruction: Store at 4°C on dry atmosphere. After reconstitution with 1 mL of deionized water, store at -20°C.

Aliquot to avoid repeated freezing and thawing.

References:

1. Protein biotinylation. Bayer EA, Wilchek M. Methods Enzymol. 1990;184:138-60.