

Datasheet

Rabbit Anti-Rat IgG F(ab')₂ secondary antibody (Peroxidase)

Catalog Number: PAB10876

Regulatory Status: For research use only (RUO)

Product Description: Rabbit anti-rat IgG recognizes rat IgG F(ab')₂ fragment. This secondary antibody was purified using antigen affinity chromatography. The antibody is conjugated with Peroxidase.

Immunogen: Rat IgG F(ab')₂ fragment

Host: Rabbit

Reactivity: Rat

Applications: Dot, EIA, ELISA, IHC, WB
(See our web site product page for detailed applications information)

Protocols: See our web site at
<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Form: Lyophilized

Conjugation: HRP

Purification: This product was prepared from monospecific antiserum by immunoaffinity chromatography, followed by solid phase adsorption(s) to remove any unwanted reactivities.

Conjugation Note:

Recommend Usage: The optimal working dilution should be determined by the end user.

Storage Buffer: Lyophilized from 5 mM sodium phosphate, 70 mM NaCl, 50% (v/v) glycerol, pH 7.4 (0.01% (w/v) gentamicin sulfate (Do NOT add sodium azide!), 10 mg/mL BSA (Immunoglobulin and Protease free))

Storage Instruction: Store at 4°C on dry atmosphere. After reconstitution with 1 mL of deionized water, store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

References:

1. Immunohistochemistry with enzyme labeled antibodies: a brief review. Farr AG, Nakane PK. J Immunol Methods. 1981;47(2):129-44.