

Datasheet

Rabbit Anti-Sheep IgG (H&L) secondary antibody (Alkaline Phosphatase)

Catalog Number: PAB10901

Regulatory Status: For research use only (RUO)

Product Description: Rabbit anti-sheep IgG recognizes sheep IgG whole molecule. This secondary antibody was purified using antigen affinity chromatography. The antibody is conjugated with Alkaline Phosphatase.

Immunogen: Sheep IgG whole molecule

Host: Rabbit

Reactivity: Sheep

Applications: Dot, EIA, ELISA, IHC, WB
(See our web site product page for detailed applications information)

Protocols: See our web site at
<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Form: Liquid

Conjugation: Alkaline Phosphatase (Calf Intestine)

Purification: This product was prepared from monospecific antiserum by immunoaffinity chromatography, followed by solid phase adsorption(s) to remove any unwanted reactivities.

Recommend Usage: The optimal working dilution should be determined by the end user.

Storage Buffer: In 50 mM Tris-Cl, 150 mM NaCl, 1 mM MgCl, 0.1 mM ZnCl₂, 50% (v/v) glycerol, pH 8.0 (0.1% sodium azide)

Storage Instruction: Store at 4°C. Do not Freeze.

References:

1. Peroxidase labelled antibody and Fab conjugates with enhanced intracellular penetration. Avrameas S, Ternynck T. Immunochemistry. 1971 Dec;8(12):1175-9.