

Datasheet

Goat Anti-Rat IgG F(ab')₂ secondary antibody (Biotin)

Catalog Number: PAB11079

Regulatory Status: For research use only (RUO)

Product Description: Goat anti-rat IgG recognizes rat IgG F(ab')₂ fragment. This secondary antibody was purified using antigen affinity chromatography. The antibody is conjugated with Biotin.

Immunogen: Rat IgG F(ab')₂ fragment

Host: Goat

Reactivity: Rat

Applications: ELISA, IHC, WB

(See our web site product page for detailed applications information)

Protocols: See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Specificity: No reaction was observed against anti-Pepsin, anti-goat IgG F(c), rat IgG F(c) or bovine, horse and human serum proteins.

Form: Lyophilized

Conjugation: Biotin

Purification: This product was prepared from monospecific antiserum by immunoaffinity chromatography, followed by solid phase adsorption(s) to remove any unwanted reactivities, pepsin digestion and chromatographic separation.

Conjugation Note:

Recommend Usage: The optimal working dilution should be determined by the end user.

Storage Buffer: Lyophilized from 5 mM potassium phosphate, 125 mM NaCl, 50% (v/v) glycerol, pH 7.6 (0.05% (w/v) sodium azide, 15 mg/mL BSA (Immunoglobulin and Protease free))

Storage Instruction: Store at 4°C on dry atmosphere.

After reconstitution with 0.1 mL of deionized water, store at -20°C or lower.

Aliquot to avoid repeated freezing and thawing.