

Datasheet

Mki67 polyclonal antibody

Catalog Number: PAB12127

Regulatory Status: For research use only (RUO)

Product Description: Rabbit polyclonal antibody raised against synthetic peptide of Mki67.

Immunogen: A synthetic peptide corresponding to amino acids 1600-1700 of mouse Mki67.

Host: Rabbit

Reactivity: Human, Mouse, Rat

Applications: Flow Cyt, IHC, IHC-P
(See our web site product page for detailed applications information)

Protocols: See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Form: Liquid

Recommend Usage: Flow Cytometry (1:100)
Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) (1:100-1:500)
The optimal working dilution should be determined by the end user.

Storage Buffer: In PBS (0.05% sodium azide)

Storage Instruction: Store at 4°C. For long term storage store at -20°C.
Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 17345

Gene Symbol: Mki67

Gene Alias: D630048A14Rik, Ki-67, Ki67