

## Datasheet

### Mki67 polyclonal antibody

**Catalog Number:** PAB12128

**Regulatory Status:** For research use only (RUO)

**Product Description:** Rabbit polyclonal antibody raised against synthetic peptide of Mki67.

**Immunogen:** A synthetic peptide corresponding to amino acids 1700-1800 of mouse Mki67.

**Host:** Rabbit

**Reactivity:** Human, Mouse

**Applications:** ICC, IF, IHC-P  
(See our web site product page for detailed applications information)

**Protocols:** See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

**Form:** Liquid

**Recommend Usage:** Immunocytochemistry (1:50-1:200)  
Immunofluorescence (1:50-1:200)  
Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) (1:100)  
The optimal working dilution should be determined by the end user.

**Storage Buffer:** In Tris-citrate/phosphate buffer, pH 7.0-8.0 (0.09% sodium azide)

**Storage Instruction:** Store at 4°C. Do not freeze.

**Entrez GeneID:** 17345

**Gene Symbol:** Mki67

**Gene Alias:** D630048A14Rik, Ki-67, Ki67